



C30-V

Your Companion, We Care

Portable & Durable for
Veterinary Monitoring



COMEN



• Intuitive Touchscreen Operation

• Transportation handle

• With only 1.5kg weight for transportation

• Command softkey board for convenient operation

Portable & Durable for Veterinary Monitoring

- 4.3-inch high resolution touchscreen display up to 2 waveforms.
- Lightweight, ergonomic design with built-in transportation handle.
- Standard features includes 3/5 lead ECG, Dual temps, CO2, SpO2, NIBP, Respiration.
- Extensive data storage capability for alarms, events, NIBP measurement and up to 48 hours trend storage & 25 minutes review.



Suntech Veterinary NIBP Cuff



Nellcor Veterinary SpO₂ Sensor



Shenzhen Comen Medical Instruments Co., Ltd.

Add: No.2 of FIYTA Timepiece Building, Nanhuan Avenue, Matian Sub-district, Guangming District, Shenzhen, 518106, Guangdong, China.
 Tel:+86-755-2640 8879 400-700-9488 Fax:+86-755-2643 1232 Website: en.comen.com E-mail: info@szcomen.com