


C50-V & C80-V

Your Companion, We Care

Portable & Durable for
Veterinary Monitoring



COMEN

- Intelligent Alarm: Auto-identification of alarm level

- Intuitive touch screen design

- Fan-less cooling design could avoid secondary air pollution and noise

- Greater data storage with NIBP trend review and alarm events

- Excellent battery capacity and performance in supporting up to 6/5.5 hours

- Support command board Keys and trim knob



C50-V/C80-V Veterinary Monitor

- Standard 10.4 inch/12 inch touch screen, convenient for operating
- Support features includes 3/5 lead ECG, Suntech NIBP Dual-IBP Dual-Temperature, RESP, Nellcor/Masimo SpO₂
- Up to 8 waveforms can displayed on the screen
- Support power-off storage and with large data storage capability for 200 alarm review, 2000 NIBP measurement and 160 hours trend storage
- Special ECG algorithms for veterinary use, accurate, reliable and able to meet the parameter measurement of animals in different body types.



Suntech Veterinary NIBP Cuff



Masimo Veterinary SpO₂ Sensor



Nellcor Veterinary SpO₂ Sensor



Veterinary 3-lead ECG Cable (Monitor)



Veterinary EtCO₂ Sampling Lines



Veterinary IBP Connector



Shenzhen Comen Medical Instruments Co., Ltd.

Add: No.2 of FIYTA Timepiece Building, Nanhuan Avenue, Matian Sub-district, Guangming District, Shenzhen, 518106, Guangdong, China.
Tel:+86-755-2640 8879 400-700-9488 Fax:+86-755-2643 1232 Website: en.comen.com E-mail: info@szcomen.com