

NEO V

Animal Dental System

Polishing, Cutting, Scaling and
Cleaning, ALL-IN-ONE



Oil Free Air Compressor

- Power: AC 120/220V 50/60Hz
- Power: 550W
- Current: 2.4A
- Airflow: 70L/min
- Exhaust Pressure : 0.8Mpa
- Noise: ≤60dB
- Temperature of usage 5-40°C
- Atmospheric Pressure: 70-106kPa



Ultrasonic Scaler Handpiece

- LED self-illumination, make the operation easier

Disinfectant system can be upgraded



**High/low speed & 3 way syringe hand pieces.
Clean /waste water bottles.
Specific suction handpiece with tube.
A magnetic rubber tray for instruments(exclude instrument).**

High-Speed Handpiece for Cutting

- Single spray
- With LED light
- Rotate speed 320000-400000rpm



Low-Speed Handpiece for Polishing

- Inner Water Spray
- Working pressure 0.3Mpa-0.35Mpa
- Rotate speed: 20000-40000rpm



3-way Syringe

- available for water,air, and spray




Detachable LED Scaler for scaling / cleaning

- Output Power: 3W-20W
- Water Pressure: 0.01Mpa-0.5Mpa



Suction Handpiece

- Take saliva, wastewater away from the pet's mouth



Low noise ≤ 60 dB, smaller tremble, smoother feel and more polishing surface



With spring-type structure and non-flying tip technology to compress the shaft, the more it turns, the tighter it becomes



Top quality bearing, ultra-long life time to handpieces , 4 holes turbine to ensure stable rotation